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The Guidelines for Using the Sterilization Room in the Department of Life Sciences

1. The sterilization autoclave in this sterilization room is exclusively available for use by all faculty and students of this department.
2. Users must comply with the usage guidelines and supervision regulations set by the department's management staff when operating equipment or engaging in activities within the sterilization room. If improper operation or behavior is observed, and after a warning the issue persists, the department reserves the right to revoke the user's permission to use the sterilization room.
3. The sterilization autoclave is a high-temperature, high-pressure device. Users must complete the department's annual laboratory hygiene and safety course and obtain approval from their supervising professor before use.
4. Undergraduate students who need to use the sterilization room for projects or science fairs must also comply with the third guideline above. They must be assisted by their supervising professor or senior colleagues in the same lab when operating the autoclave. Users who fail to follow these regulations and cause damage to property or personnel will be held responsible.
5. Before using the autoclave, verify that it is not already in use. Fill out the operation log before and after use. If any abnormalities are found, immediately inform the laboratory manager.
6. Ensure that the pressure and temperature gauges are within the maximum specified limits according to the user manual. Always operate the autoclave within the normal range. Before sterilization, remember to add water (to maintain the cleanliness of the autoclave, do not use tap water).
7. After the sterilization process is complete, ensure that the pressure has returned to zero and the temperature has dropped before opening the autoclave door. Do not forcibly open the autoclave if the pressure has not reached zero or the temperature is still above 65°C, to avoid accidents.
8. Waste materials should be sterilized separately from other items. This sterilization room is only available for the sterilization of waste materials in autoclaves numbered 5-6.
9. To prevent breakage of items or leakage of liquids during sterilization, do not overload the autoclave. Each load should not exceed 80% of the total capacity of the autoclave. Ensure that the top of the items is at least 3 inches (about 8 cm) away from the autoclave wall, and leave gaps between each package. Do not pack items tightly together. Items should not be placed directly on the bottom of the autoclave; they must be placed in stainless steel containers with handles before being loaded.
10. After use, users must clean both the inside and outside of the autoclave. All sterilized items should be removed the same day; do not leave them in the sterilization room.
11. Any violation of the above regulations will be reported by the sterilization room manager to the supervising professor. The violating user's laboratory will be suspended from using the sterilization room for one week. During the suspension, the laboratory is responsible for cleaning the sterilization room. If the same laboratory commits further violations, the suspension will be extended to two weeks, and a fine of NT\$5,000 will be imposed (deducted from the professor's research fund). If a laboratory violates the regulations three or more times in the same academic year, the case will be submitted to the department's academic affairs meeting for a decision.